


**How to Use the
GILBERT MICROSCOPE**

**Cover
Table of Contents
pp 1-19**

HOW TO USE THE

Gilbert

MICROSCOPE



ANOTHER GILBERT HALL OF SCIENCE PRODUCT

EXPLORING THE WORLD WITH THE MICROSCOPE

Edited by

OSCAR W. RICHARDS, Ph.D.

University of Oregon, B.A. 1923—M.A. 1925

Yale University, Ph.D. 1931

In collaboration with

ALFRED C. GILBERT, M.D.

Yale University 1909



Copyright 1936

The A. C. Gilbert Company

Erector Square, New Haven, Connecticut, U.S.A.

DEVELOPED AT THE GILBERT HALL OF SCIENCE

M2987

Printed in U.S.A.

LIST OF EXPERIMENTS THAT CAN BE DONE WITH THE No. 8 SET

1	2	3	4	5	6	7	8	9	10	11	12
13	14	22	23	31	32	33	34	35	36	37	38
39	40	41	42	43	44	45	46	47	48	49	50
51	53	54	55	56	57	58	59	64	66	67	68
79	80	81	82	87	88	89	92	94	95	96	97
98	100	101	102	109	110	112	115	116	117	118	119
120	121	122	123	124	126	127	128	129	130	133	135
136	137	138	140	141	142	143	144	145	146	147	148
149	150	151	153	154	155	156	157	158	160	161	165
166	167	168	169	172	173	174	175	176	177	178	179
180	181	182	183	184	185	186	187	189	192	194	195
196	197	200	202	203	204	205	206	207	208	209	210
211	212	214	215	216	217	218	219	220	221	222	223
224	225	226	227	228	229	230	231	233	234	235	236
238	239	241	243	245	246	247	248	249	250	251	252
253	254	255	256	257	258	259	260	263	265		

TABLE OF CONTENTS

Introduction	Page
Warning	14
Instructions for cleaning and correctly replacing lenses in your Microscope....	15
	16

CHAPTER 1

THE SIMPLE MICROSCOPE

EXPERIMENT

1. Using the Hand Lens	17
2. Magnifying Power	17
3. Magnifying Power of a Lens	17
4. Making a Leuwenhoek Microscope	18
5. Cloth	18
6. Finger Tip Patterns	18
7. Finger Prints	18
8. Lace	18
9. How Printing Appears	18
10. Printed Pictures	19
11. Some Interesting Experiments	19
12. How to use the Dissecting Microscope	19

CHAPTER 2

THE COMPOUND MICROSCOPE

13. Reversal of Image by the Compound Microscope	24
14. Making a Temporary Slide	25
15. The Number One Objective	25
16. The Number Two Objective	25
17. The Number Three Objective	26
18. Adjusting the Mirror	26
19. The Draw Tube	26
20. Magnification	26
21. The Importance of the Cover Glass	27
22. Locating Dirt on the Microscope Lenses	27

CHAPTER 3

CARE OF THE EYES WHEN USING THE MICROSCOPE

23. Determining Whether or not your Glasses are to Correct for Astigmatism	28
--	----

CHAPTER 4

OBSERVATIONS AND RECORDS

CHAPTER 5 COLLECTING AND STORING MATERIAL

CHAPTER 6 MEASURING WITH THE MICROSCOPE

	Page
EXPERIMENT	
24. The size of the field seen with the Microscope	31
25. Measuring a Salt Crystal	32
26. Vertical Measurement with the Microscope	32
27. Thickness of a Paint Film	32
28. Estimating the Enlargement of a drawing	33

CHAPTER 7 PHOTOMICROGRAPHY

29. Making Microphotographs	34
---------------------------------------	----

CHAPTER 8 HAIR AND THE FUR INDUSTRY

30. A Brunette Human Hair	35
31. A Blond Hair	35
32. A Gray Hair	35
33. A Red Hair	35
34. Straight Hair	35
35. Curly Hair	36
36. Permanent Waving of Hair	36
37. Examine a Bit of Wool	36
38. A Cat Hair	36
39. A Dog Hair	37
40. A Rabbit Hair	37
41. Identifying a Hair	37

CHAPTER 9 MAKING PERMANENT SLIDES

42. Making a Permanent Slide with Karo Corn Syrup	39
43. Making a Permanent Slide with Balsam	38

CHAPTER 10 HOW WE SEE WITH THE MICROSCOPE

44. Colored Image	39
45. Refraction Images	39

EXPERIMENT

46. Determining the Direction of a Refraction Image	40
47. The Becke Line	40
48. Determining Refractive Index	40
49. Refractive Index of a Heavy Oil	41
50. Refractive Index of Clove Oil	41
51. Refractive Index of Oil Bubbles	41
52. Opaque Materials	41
53. Color from Diffraction Patterns	42

CHAPTER 11 PREPARATION OF OBJECTS

54. Sifting	43
55. Crushing and Sifting	43
56. Kneading	43
57. Sedimentation	43
58. Separation by Special Mixtures	44
59. How to use the Collecting Tube	44
60. Obtaining Rice Facing	44
61. Starch or Tale?	44
62. Sampling Soil	44
63. Clearing with Oil	45
64. Mounting in Oil	45
65. Filtration	45
66. Teasing Paper Specimens	46
67. Teasing Alent Fibers	46
68. Teasing Potato Cells	46
69. Sectioning	46
70. Use of the Microtome	47
71. Cutting Thin Materials	47
72. Sectioning with the Microtome	47

CHAPTER 12 FUNDAMENTAL UNITS OF THE INORGANIC WORLD

73. Dendrite Crystals of Aspirin	48
74. Starch Test with Iodine	48
75. Crystal Growth	48
76. Crystal Growth from a Metal	49
77. Supersaturation	49
78. Speeding Crystal Growths	49
79. Water of Crystallization	49
80. Water in Washing Soda Crystals	49
81. Broken Crystals	50
82. Crystal Repair	50
83. Precipitation Membrane Formation	50
84. Crystals from Fusion	50
85. Effect of Gelatine on Crystal Growth	50
86. The Effect of Albumen on Crystal Growth	50

CHAPTER 13 THE CLASSIFICATION OF CRYSTALS

EXPERIMENT		Page
87. Cubic Crystals of Salt	51	51
88. Octahedral Crystals of Alum	51	51
89. Tetragonal Crystals	52	52
90. Rhombic Crystals of Potassium Nitrate	52	52
91. Monoclinic Crystals of Borax	52	52
92. Triclinic Crystals of Boric Acid	52	52

CHAPTER 14 CHEMICAL MICROSCOPY

93. Solution in Other Liquids than Water	53
94. Separation of Salt by Solution and Filtration	53
95. Test for Salts from Solution	54
96. Stirp Method of Filtering	54
97. Testing Tooth Paste	54
98. Salts in Soil	54
99. Another Test for Tooth Paste	54
100. Micro-Decantation	54
101. Sublimation of Iodine	55
102. Sublimation of Tannin from Tea	55
103. Testing a Raisin for Copper	55
104. Testing Apple Seeds for Iron	55
105. Testing Lettuce for Iron	55
106. Testing Ash of Bean for Iron	56
107. Testing Peas for Copper	56
108. Reagents on a Fiber	56
109. Test for a Carbonate	56
110. Testing Clam Shell for Carbonate	56
111. Hydrogen Gas from Zinc	56
112. Carbonate from Expired Air	56
113. Proving a Gas to be Carbon Dioxide with Lime Water	57

CHAPTER 15 FUNDAMENTAL UNITS OF LIVING ORGANISMS

114. Cork Cells	58
115. Cells from Celery	58
116. Plasmolysis of a Cell	58
117. Sugar and Salt Solutions	58
118. Chloroplasts	59
119. Animal Cells	59
120. The Form of Soap Bubbles	59
121. The Cells of Cucumber	60
122. The Cells in Onion	60
123. Oil Stored in Cells	60
124. Staining Cells	61

EXPERIMENT		Page
125. Staining with Saffranin	61	61
126. Human Red Blood Cells	61	61
127. Red Blood Cells of a Frog	61	61
128. Cells from Lean Meat	61	61
129. The Fibers in a Nerve	61	61
130. Plant and Animal Cells	61	61

CHAPTER 16 UNICELLULAR PLANTS AND ANIMALS

131. The Paramecium	65
132. Speed of a Paramecium	65
133. Defences of a Paramecium	65
134. The Euglena	65
135. Comparison of the Euglena with the Paramecium	65
136. The Amoeba	66
137. What Chalk is made of	66
138. Foraminifera	66
139. Diatoms	66
140. Diatom Skeletons	66
141. Testing Silver Polish	68
142. Testing Tooth Paste Abrasive	68
143. Living Diatoms	68
144. Spirogyra	68
145. Oscillatoria	68
146. Pond Scum	68
147. Volvox	68
148. Round Bacteria	69
149. Rod-Shaped Bacteria	69
150. Spiral Bacteria	69
151. Bacteria Slide with Thionine	70
152. Bacteria Slide with Saffranin	70
153. Yeast Cells from a Yeast Cake	70
154. Starch in a Compressed Yeast Cake	70
155. Growing Yeast Plants	70
156. Bread Mold	71
157. Mucor	71
158. Molds of Fruit and Jelly	71

CHAPTER 17 MOSESSES, FERNS AND FLOWERING PLANTS

159. Moss Plants	73
160. Ferns	73
161. Examination of a Root from a Seed	73
162. Sections of Roots	73
163. Wood Fibers	76
164. The Structure of a Grass Stem	76

EXPERIMENT		Page
165. A Leaf	76	76
166. The Structure of the Flower	77	77
167. Pollen Grains	77	77
168. Pollen Tube Growth	77	77
169. Starch from Flour	77	77
170. The Effect of Cooking in Starch	77	77
171. The Digestion of Starch	78	78

CHAPTER 18 THE MANY CELLED ANIMALS

172. Hydra	79
173. A Round Worm	79
174. A Flat Worm	79
175. Regeneration of an Injured Worm	79
176. Crustaceans	82
177. Sea Moss	82
178. Wheel Animals	83
179. The House Fly	83
180. Mosquito	83
181. Fruit Fly	83
182. The Cause of the Cricket's Chirp	83
183. The Honey Bee	83
184. Insect Wings	84
185. Scales of a Butterfly's Wings	84
186. Giliated Cells from a Clam	84
187. Fish Scales	85
188. Why Feathers are Stiff	85
189. The Development of a Frog Egg	87

CHAPTER 19 MAKING PERMANENT MOUNTS OF MOIST TISSUES

190. Prepare Some Slides of Stained Tissues	88
---	----

CHAPTER 20 PRINTING

191. Engraving	89
192. A Line Cut Impression	89
193. A Coarse Screen Halftone Picture	89
194. A Fine Screen Picture	89
195. Printed Color Plate	90
196. Rotogravure Pictures	90
197. A Benday Picture	90

CHAPTER 21 THE STRUCTURE OF THE PHOTOGRAPHIC IMAGE

EXPERIMENT		Page
198. Structure of a Photographic Negative	91	91
199. Limitations of Enlargement	91	91
200. Correspondence in Grain Between Negative and Enlarged Print	91	91
201. Fine Grain Developer	91	91
202. Reversal Film	92	92
203. Natural Color Films	92	92

CHAPTER 22 METALLURGY AND THE MICROSCOPE

204. Rusty Iron	93
205. Metal Surfaces	93
206. Metal Plated Surfaces	93
207. Flaws in a Broken Surface	93
208. Alloys and Its Components	93

CHAPTER 23 TEXTILE MICROSCOPY

209. Wool	95
210. Shoddy	95
211. Mohair	95
212. Cotton	95
213. Part Wool Cloth	95
214. Mercerized Cotton	95
215. Linen	96
216. Natural Silk	96
217. Colloidion Silk	96
218. Cuprammonium Silk	97
219. Viscose Silk	97
220. Acetate Silk	97
221. Sheet Cellulose Products	97
222. Rope Fibers	97

CHAPTER 24 THE USE OF THE MICROSCOPE IN THE PAPER INDUSTRY

223. Chemical and Mechanical Paper Pulp	98
224. Coniferous and Non-Coniferous Paper Pulp	98
225. Unbleached Paper	99
226. Cardboard and Strawboard	99
227. Linen Paper	99
228. Paper Sizing	99

CHAPTER 25 OTHER APPLICATIONS OF THE MICROSCOPE

EXPERIMENT	Page
229. Sound Groove in a Phonograph Record	100
230. Spices and their Adulteration	100
231. Parts of Insects	100
232. Designs of Diatoms	100
233. Cigar and Cigarette Ash Collections	100
234. Snow Flakes	100

CHAPTER 26 DETECTIVE WORK WITH THE MICROSCOPE

235. Identifying Soil Samples	102
236. The Vacuum Cleaner Detective	102
237. The Handwriting Expert	102
238. Typewriter Detection	103
239. Can Hand Writing be Disguised?	103
240. Safety Paper	103
241. Oil Paintings	104
242. Kinds of Finger Prints	104
243. The Home Detective	104
244. Measuring the Magnifying Power of the Compound Microscope	106

CHAPTER 27 THE USE OF POLARIZED LIGHT

245. Polarization by Reflection	108
246. The use of two Polarizing Filters	108
247. Using Polarized Light with the Microscope	109
248. Potato Starch in Polarized Light	109
249. Starches and Polarized Light	109
250. Identifying Starch with Polarized Light	109
251. Plant and Animal Hairs with Polarized Light	109
252. Polarization Images from Wood	110
253. Sections of Shells, Horn, Bone and Teeth	110
254. Fish Scales	110
255. Cotton, Silk, Wool and the Natural Fibers	110
256. Separating the Artificial Silks with Polarized Light	110
257. Cellulosic Wrapping Materials	110
258. Retardation Tint Plates	110
259. Effect of Rotating the Analyzer	111
260. Examination of Minerals with Polarized Light	111
261. Doubly Refractive Colors of some Crystals	111
262. Crystals from Fusion	111
263. Crystals from Sublimation	111
264. Converting the Microscope into a Polarimeter	111
265. Inverting Sugar	112
266. Making Zinc Acetate Crystals	112

Bibliography	Page
Books on the Microscope and its use	113
Books devoted to methods	113
Books on technology	113
Books on biology	113
List of Replacement Parts Purchaseable at The A. C. Gilbert Co.	114
Sources of Supplies not handled by The A. C. Gilbert Co.	115

INTRODUCTION

Now that you own a Microscope you are equipped to undertake a journey of exploration such as you probably have never before enjoyed, for touring the world with your microscope for hidden sights is one of the most exciting pleasures that you can imagine.

Just think of the thrills there are in store for you.

Hidden mysteries and secrets of nature and science are brought to your clear vision through powerful, searching lenses and one after another amazing and unbelievable secrets are revealed.

You will see hundreds of new things — too small to find with just your eyes and you will marvel at intricate details of how things are made and put together.

You'll not only explore new and fascinating worlds but learn how to make tests, know about plants, animals, fur, textiles, paper, crystals, and do detective work as well.

The microscope, especially in this modern age is recognized as indispensable and few industries can get along without its aid.

Physicians and their assistants use it for blood counts and in examining bacteria which causes disease, scientists in making new and startling discoveries and you too can find no end of interesting objects by exploring your home and the things you and your folks use each day.

Each season of the year offers exciting experiments of its own and many of the experiments described in this manual may be performed anytime, even in mid-winter.

You will find some things which show nicely under the microscope with little care, but to really use the microscope you will wish to know how it works and how to prepare objects for examination. Even if you have used a microscope before, you will enjoy reading the first chapters of the manual and remember to make discoveries properly you should follow the directions given. Very few objects can be looked at with the microscope without some special treatment and this manual is written to tell you how to have the most fun with your new equipment. No manual, however long, could tell you of all the interesting things to see, but you will find thousands of objects all about you, at all times of the year, wherever you may be. Many people find that the microscope becomes a pleasant companion and make microscopy a delightful hobby. Years ago men either had to make their microscopes or had to pay a great deal for them, but now any boy may have a good microscope and can explore with it. Some of you may find that the microscope is necessary to your life work

WARNING!

This set is not intended for children who cannot read and understand the accompanying Instruction Book.


This set does not contain dangerous poisons and the chemicals mentioned in this manual are not embraced under the term "poisons." They are perfectly safe to use if handled carefully and intelligently. They are not intended to be taken by mouth or swallowed, and no intelligent person would be expected to use them for such purposes. It is necessary, however, to emphasize the fact that carelessness on the part of the experimenter can always lead to trouble. The author suggests, therefore, that all experimentation be carried out cautiously and according to the directions, especially when manipulations like heating is involved, or when gases are evolved in the reactions.

Before performing any experiments outlined in the manual, the following instructions should be read carefully and observed.

Before performing experiments, be sure to spread a thick layer of newspapers or other protective material over the table, so that hot liquids, candle grease, etc., will not injure the table.

Always read an experiment entirely through before starting to perform it. By following this rule many mistakes may be avoided.

Never point the open end of a test tube, while heating, at yourself or anyone nearby, as it may suddenly boil over, causing burns or injuring clothing. For the same reason never smell at the open end of a test tube while heating, or put your face near it.



INSTRUCTIONS FOR CLEANING AND CORRECTLY REPLACING LENSES IN YOUR MICROSCOPE

Microscope S6 has a collective lens and diaphragm located in main barrel and spaced as shown in diagram. After cleaning lenses it is important that spacing tubes, etc. be replaced exactly as shown.

Microscope S15 lenses are mounted in a short metal tube which can be taken from the large tube. Remove the eye lens in its mount, push long spacing tube up with finger and remove short spacing tubes and lens as you push them up. Clean lenses and replace in the order shown in diagram. For further information on the care of microscope lenses see Experiment 22 on Page 27, entitled, "Locating Dirt on the Microscope Lenses."

The Simple Microscope

Any round piece of transparent material may be used to make a simple microscope, or hand magnifying glass. No one knows who discovered the first one, but we do know that many years ago clear pebbles were used to make objects look larger. The clearer the stone the better one can see through it. Glass is used now to make magnifying glasses. The glass may be round, or one side may be flat and the other side round. All that is necessary is that the center of the glass be thicker than the edge.

EXPERIMENT 1. Using the Hand Lens

Pull the simple lens from the support rod of the dissecting microscope of your set, fig. 1, and hold it near a piece of newspaper or other object. Bring the lens toward your eye until you find the position that gives the best view. Are the letters, or is the object right side up? Move the object while you look at it with the lens and note whether it appears to move in the same direction.

EXPERIMENT 2. Magnifying Power

Look at a hair or a small bit of wire. Compare the size of it as seen through the lens with the actual size. About how many times larger do you think it makes the hair?

EXPERIMENT 3. Magnifying Power of a Lens

A more accurate method of measuring the magnifying power is to determine the focal length of the lens. Hold a piece of white paper so that it faces a window or a light. Place the lens in front of the paper and move it slowly toward the window or light until you see a picture of the window or light on the paper. Then measure the distance from the center of the lens to the paper and this distance is the focal length. Divide 10 by the number of inches focal length and the result is the magnifying power. The lens that I have, has a focal length of 1 and $\frac{1}{2}$ inches. Therefore its magnifying power is $10 \div 1.5$ which is just a little over six and one-half times. Everything that I look at with it appears to be nearly seven times larger.

The magnifying power depends on how curved the lens is. The greater the curvature the larger appears the object seen through the lens. Smaller lenses have greater curvature. Tiny glass beads or round drops of water make high powered lenses. Leeuwenhoek (pronounced Lay'wen hook) made many simple microscopes by grinding very small glass beads and mounting them in a hole between two pieces of metal (fig. 2). Because it was hard to hold such a microscope steady he added an arm to hold the object and arranged screws to move the object and to hold it in place. Then the microscope could be passed around and his friends could see the interesting things. He saw bacteria, small plants and animals and his discoveries were




Fig. 1




Fig. 2

from a model made by the Bausch and Lomb Optical Co. in 1932 to celebrate the three hundredth anniversary of his birth.

EXPERIMENT 4. Making a Leeuwenhoek Microscope

Take a piece of cardboard and make a small hole through it with a nail. The hole should be less than $\frac{1}{8}$ inch in diameter. If the cardboard is soft rub a little warm candle grease or paraffine into it to make it water proof. With a medicine dropper place a small drop of water in the hole. This is the lens. The drop should be small enough to round up. Place a fly on a pin and hold it close to the drop and look at it through the drop. You will find that the fly has to be held closer to the drop than it did to the hand lens. About how much does this simple microscope magnify? When the water evaporates it will have to be replaced. By using different sized holes and drops of water you can make similar microscopes with different powers.

The simple microscope, or hand lens, is used for preliminary examination of objects. The following experiments describe some of the observations you may make and you will think of many others.

EXPERIMENT 5. Cloth

Examine some cloth, such as a handkerchief and see how it was woven. How many threads are used to one-quarter inch? Do cloths which feel finer have more or less threads to the quarter inch?

EXPERIMENT 6. Finger Tip Patterns

Examine the ridges on your finger tips with the hand lens. Do they all have the same pattern? Are they the same on your friend's fingers?

EXPERIMENT 7. Finger Prints

If you have a stamp pad place a finger on it to cover it with ink and then place it on a piece of paper to give a finger print. If you do not have a stamp pad rub off some lead from a pencil and use it instead of the pad. Compare different finger prints using your lens.

EXPERIMENT 8. Lace

Look at some lace and see if it is woven the same as cloth. To see the object clearly it must be held so that it is in the light and not shaded by the lens. Take your object to a window or to a light and you will soon learn how to hold

the lens and object for best vision. Sometimes one can see better when the lens is held close to the eye.

EXPERIMENT 9. How Printing Appears

Look at different kinds of printing to see if the letters are equally clear. Look at a bit of colored "funny" paper and see how the colors are printed. Some will be plain color and others made by printing dots of a different color over the first color. The dots seem to blend with the plain color when seen by the unaided eye.

EXPERIMENT 10. Printed Pictures

Compare pictures made on good paper in a magazine with those printed in a newspaper. You will see that the picture is actually printed as dots and the coarser the paper, like newspapers, the coarser are the dots and the less clear the picture details become.

EXPERIMENT 11. Some Interesting Experiments

Minerals, frost on the window pane, crystals, dust, etc. will be interesting subjects when seen under your hand lens. In fact any object that can be brought into the light becomes more interesting under the magnifying glass. Be sure to look at plants, flowers, insects, etc.

EXPERIMENT 12. How to Use the Dissecting Microscope

If you wish to look at part of an insect it is convenient to have a stand and a rod for holding the lens so that your hands are free. For this work one puts the lens back onto the support rod of the dissecting microscope (fig. 1). The insect is put on one of the glass slides placed under the lens. The lens is slid up or down on the rod until the insect is seen clearly, or as we say, the lens is focused on the insect. The rod will hold it so you can move the insect around with the dissecting needles (fig. 3B) or take it apart with the forceps (sometimes called tweezers, fig. 3C). This dissection microscope will be used to prepare material for examination with the compound microscope as will be described later.




Fig. 3